



1743 Creekside Drive, Suite 200  
Folsom, CA 95630

530.626.7000 | [www.cgt.global](http://www.cgt.global)

email: [info@cgt.global](mailto:info@cgt.global)

## Cryopreserved Product

Mobilized Peripheral Blood  
CD34+ Stem/Progenitor Cells, G-CSF

|          |              |                   |
|----------|--------------|-------------------|
| Catalog# | MLEG34000.2C | 0.2 million cells |
|          | MLEG34001C   | 1 million cells   |
|          | MLEG34005C   | 5 million cells   |
|          | MLEG34010C   | 10 million cells  |
|          | MLEG34020C   | 20 million cells  |

## PRODUCT DESCRIPTION

Human Mobilized Peripheral Blood CD34+ Stem/Progenitor Cells are positively selected from mobilized peripheral blood mononuclear cells.

Mobilized Leukopaks<sup>®</sup> are collected using the Spectra Optia<sup>®</sup> Apheresis System from donors that are injected with 10µg/kg/day Granulocyte-Colony Stimulating Factor (G-CSF). Mononuclear cells are enriched from the Mobilized Leukopak using an RBC lysis protocol. CD34+ hematopoietic stem cells (HSCs) are then selected for using immunomagnetic anti-CD34 microbeads from the mononuclear cell pool, leaving a highly pure population of CD34+ stem and progenitor cells.

Cells were obtained using Institutional Review Board (IRB) approved consent forms and protocols.

## CRYOPRESERVATION

Cryopreserved products allow for prolonged storage before use. Cell products contain 10% DMSO to minimize cell death during freezing. All cryopreserved products are stored in containers designed and tested for ultra-low temperatures at long time intervals. We normally ship cryopreserved items on dry ice, but can also use a cryoshipper at the customer's request.

## SAMPLE COLLECTION AND PROCESSING

All samples are collected on-site at our Stem Cell Collection Center, or at nearby partner hospitals or clinics. Apheresis donors are transfused with ACD-A during the collection process. Samples are then quickly processed in our on-site laboratory to achieve maximum viability and quality. Cryopreserved cells are frozen at -1°C/minute in a -80°C freezer, and then transferred to liquid nitrogen.

Infectious disease testing for HIV, HBV, and HCV is performed on a sample of donor blood. Only samples with negative results within 90 days of collection are shipped unless approved by the customer. All testing is performed by a CLIA-certified lab.

## FORMAT

Stem/progenitor cells are frozen in StemSpan<sup>™</sup> + 10% DMSO. We can also use freezing media as specified by the customer.

## STORAGE

Cryopreserved cells should be maintained at -135°C or colder (in liquid nitrogen). The cells are warranted for twelve months from the date of receipt if stored at -135°C or colder. Storage of cells at -80°C for less than one month should maintain cell viability but is not covered by the warranty.

## THAWING INSTRUCTIONS FOR CELL PRODUCTS

Refer to the Technical Support page at [cgt.global/technical-support/](http://cgt.global/technical-support/) to access our online Thawing Protocols.

## WARNING

This product contains human tissue or other biological material and MUST be handled at Biosafety Level 2 or higher. All biological products should be treated as potentially infectious or contaminated material, even if infectious disease screening reports are negative. Follow universal precautions and wear appropriate personal protective equipment.

## PRODUCT WARRANTY

For our product warranty, please review our Terms and Conditions at [cgt.global/terms-and-conditions/](http://cgt.global/terms-and-conditions/).

FOR RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.