

StemExpress Frozen Leukopaks

Time is of the essence when it comes to your experiments. Viable source material is crucial to the success of your projects, meaning delays and logistical issues can not only be costly, but they can compromise downstream project success. StemExpress frozen Leukopaks allow you to keep the high-quality cells you need on hand – ready to thaw and use exactly when you need them. Our Leukopaks are expertly cryopreserved within 24 hours of collection to maintain cell quality and mitigate cell death that can occur in the days following the donor collection.

- ✓ Donors follow FDA 21 CFR 1271 guidelines
- Average cell viability 96.33% post-thaw as compared to 97.56% for a fresh Leukopak
- Customizable donor selection process, infectious disease screening and timepoints
- Customizable bag sizes and cell concentrations available to fit your process





Frozen Leukopaks

Cell-rich Leukopaks are collected from IRB consented donors using the Spectra Optia® Apheresis System at our FDA-registered StemExpress Collection Center. Frozen Leukopaks are cryopreserved within 24 hours of collection to ensure the highest quality product.


- ✓ 1x10e10 TNCs/full bag (prior to cryopreservation)
- √ ≥95% viability (prior to cryopreservation)
- ✓ Highly-enriched source of PBMCs


- ✓ Retention vials
- ✓ Low granulocyte/hematocrit
- ✓ HBV, HCV, HIV negative



Our frozen Leukopaks reduce time constraint barriers and improve efficiency & logistics.

- Preserves cell viability during global shipping
- Provides convenience in the timing of experiments
- Allows future access to same lot samples for additional studies
- Ideal for developing intra- and interlaboratory proficiency testing in clinical trials





Customization

We provide customizable options to ensure you receive the Leukopak you need, when you need it.

- Recallable donors
- Donor demographics & lifestyle characteristics
- Medical history
- HLA, RhD, and ABO typing available
- Additional infectious disease screening
- Collection scheduling
- Collection protocols, such as half fresh and half frozen

Leukopak Product list

Part Code	Description	Format	Size	Container
LE002.5C	Peripheral Blood 1/4 Leukopak	Cryopreserve	2.5 x 10E9	1 Bag
LE002.5F	Peripheral Blood 1/4 Leukopak	Fresh	2.5 x 10E9	1 Bag
LE005C	Peripheral Blood 1/2 Leukopak	Cryopreserve	5.0 x 10E9	1 Bag
LE005F	Peripheral Blood 1/2 Leukopak	Fresh	5.0 x 10E9	1 Bag
LE010C	Peripheral Blood Full Leukopak	Cryopreserve	1.0 x 10E10	1 Bag
LE010F	Peripheral Blood Full Leukopak	Fresh	1.0 x 10E10	1 Bag
CUSLPC01	Custom Leukopak Collection	Fresh or Cryopreserved	Custom	Custom

ORDERING & SHIPPING 530.303.3828 | 888.415.4215 sales@stemexpress.com

