



1743 Creekside Drive, Suite 200
Folsom, CA 95630

530.626.7000 | www.cgt.global

email: info@cgt.global

Fresh Product

Whole Peripheral Blood Vacutainer®
Tube Collection, Na-Heparin

| | | | | |
|----------|-----------|--------|-----------|--------|
| Catalog# | PBNAH010F | 10 mL | PBNAH110F | 110 mL |
| | PBNAH020F | 20 mL | PBNAH120F | 120 mL |
| | PBNAH030F | 30 mL | PBNAH130F | 130 mL |
| | PBNAH040F | 40 mL | PBNAH140F | 140 mL |
| | PBNAH050F | 50 mL | PBNAH150F | 150 mL |
| | PBNAH060F | 60 mL | PBNAH160F | 160 mL |
| | PBNAH070F | 70 mL | PBNAH170F | 170 mL |
| | PBNAH080F | 80 mL | PBNAH180F | 180 mL |
| | PBNAH090F | 90 mL | PBNAH190F | 190 mL |
| | PBNAH100F | 100 mL | PBNAH200F | 200 mL |
| | | | PBNAH250F | 250 mL |

PRODUCT DESCRIPTION

Human Whole Peripheral Blood is collected in 10 mL sterile vacutainer(s) containing the anticoagulant Na-Heparin.

Fresh products have a high viability without the detrimental effects of freezing, thawing, and exposure to cryoprotectants.

Blood is obtained using Institutional Review Board (IRB) approved consent forms and protocols.

SAMPLE COLLECTION AND PROCESSING

All samples are collected on-site at our Stem Cell Collection Center.

Infectious disease testing for HIV, HBV, and HCV is performed on a sample of donor blood. Only samples with negative results within 90 days of collection are shipped unless approved by the customer. All testing is performed by a CLIA-certified lab.

FORMAT

Fresh whole peripheral blood is shipped at ambient temperature.

Specific containers and media can also be prepared as requested by the customer.

STORAGE

Fresh products should be used or processed immediately upon receipt. The warranty only covers items whose specifications are tested at the time they are received.

WARNING

This product contains human tissue or other biological material and MUST be handled at Biosafety Level 2 or higher. All biological products should be treated as potentially infectious or contaminated material, even if infectious disease screening reports are negative. Follow universal precautions and wear appropriate personal protective equipment.

PRODUCT WARRANTY

For our product warranty, please review our Terms and Conditions at cgt.global/terms-and-conditions/.

FOR RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.