Stemexpress[®] 1743 Creekside Drive, Suite 200 Folsom, CA 95630

advancing science

1743 Creekside Drive, Suite 200 Folsom, CA 95630 530.626.7000 email: info@stemexpress.com www.stemexpress.com

Cryopreserved Product

Cord Blood CD133+ Stem Cells, Pooled Donors

Catalog#	CBP13300.5C	0.5 million cells
	CBP13301C	1.0 million cells
	CBP13302C	2.0 million cells
	CBP13302.5C	2.5 million cells

Product Description

Human Umbilical Cord Blood CD133+ Stem/Progenitor Cells are positively selected from umbilical cord blood mononuclear cells.

Whole umbilical cord blood is needle aspirated from the umbilical cord vein using a cord blood collection bag containing 35 mL of the anticoagulant citrate phosphate double dextrose (CP2D). Mononuclear cells are enriched from the cord blood using a density gradient centrifugation protocol. CD133+ hematopoietic stem cells (HSCs) are then selected using immunomagnetic anti-CD133 microbeads from the mononuclear cell pool, leaving a highly pure population of CD133+ stem and progenitor cells.

Single and mixed donor CD133+ stem cells are available upon request. Mixed donor pools are from two or more donors and are available in larger lots. Single donor cells are obtained from a single umbilical cord sample.

Cells were obtained using Institutional Review Board (IRB) approved consent forms and protocols.

Cryopreservation

Cryopreserved products allow for prolonged storage before use. Cell products contain 10% DMSO to minimize cell death during freezing. All cryopreserved products are stored in containers designed and tested for ultra-low temperatures at long time intervals. We normally ship cryopreserved items on dry ice, but can also use a cryoshipper at the customer's request

Sample Collection and Processing

All samples are collected at nearby partner hospitals or clinics. Samples are then quickly processed in our on-site laboratory to achieve maximum viability and quality. Cryopreserved cells are frozen at -1°C/minute in a -80°C freezer, and then transferred to liquid nitrogen.

Infectious disease testing for HIV, HBV, and HCV is performed on a sample of cord blood by a CLIA-certified lab.

Format

Isolated stem/progenitor cells are frozen in StemSpan[™] + 10% DMSO. We can also use freezing media as specified by the customer.

Storage

Cryopreserved cells should be maintained at -135°C or colder (in liquid nitrogen). The cells are warrantied for twelve months from the date of receipt if stored at -135°C or colder. Storage of cells at -80°C for less than one month should maintain cell viability but is not covered by the warranty.

StemSpan is a trademark of StemCell Technologies Inc.

Thawing Instructions for Cell Products

Refer to our "How to thaw StemExpress primary cells for optimal viability?" under our Frequently Asked Questions at stemexpress.com/faqs/ to access our online Thawing Protocol.

Warning

This product contains human tissue or other biological material and MUST be handled at Biosafety Level 2 or higher. All biological products should be treated as potentially infectious or contaminated material, even if infectious disease screening reports are negative. Follow universal precautions and wear appropriate personal protective equipment.

Product Warranty

For our product warranty, please review our Terms and Conditions at stemexpress.com/terms-and-conditions/.

FOR RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.

