

1743 Creekside Drive, Suite 200 Folsom, CA 95630

530.626.7000 | www.cgt.global

email: info@cgt.global

Cryopreserved Product

Cord Blood CD8+ Cytotoxic T Cells

Catalog# CB08005C 5 million cells

CB08010C 10 million cells
CB08020C 20 million cells
CB08025C 25 million cells

PRODUCT DESCRIPTION

Human Umbilical Cord Blood CD8+ Cytotoxic T Cells, also referred to as cytotoxic T cells (CTLs), are negatively selected from umbilical cord blood mononuclear cells .

Whole umbilical cord blood is needle aspirated from the umbilical cord vein using a cord blood collection bag containing 35 mL of the anticoagulant citrate phosphate double dextrose (CP2D). Mononuclear cells are enriched from cord blood using a density gradient centrifugation protocol. Cells expressing CD4, CD14, CD16, CD19, CD20, CD36, CD56, CD66b, CD123, $TCR\gamma/\delta$ and CD235a are depleted from the mononuclear cell population using immunomagnetic particles leaving purified, untouched CD8+ cytotoxic T cells.

Cells were obtained using Institutional Review Board (IRB) approved consent forms and protocols.

CRYOPRESERVATION

Cryopreserved products allow for prolonged storage before use. Cell products contain 10% DMSO to minimize cell death during freezing. All cryopreserved products are stored in containers designed and tested for ultra-low temperatures at long time intervals. We normally ship cryopreserved items on dry ice, but can also use a cryoshipper at the customer's request

SAMPLE COLLECTION AND PROCESSING

All samples are collected at nearby partner hospitals or clinics. Samples are then quickly processed in our on-site laboratory to achieve maximum viability and quality. Cryopreserved cells are frozen at -1°C/minute in a -80°C freezer, and then transferred to liquid nitrogen.

Infectious disease testing for HIV, HBV, and HCV is performed on a sample of cord blood by a CLIA-certified lab.

FORMAT

Isolated cells are normally frozen in CryoStor $^{\text{TM}}$ CS10 (10% DMSO). We can also use freezing media as specified by the customer.

STORAGE

Cryopreserved cells should be maintained at -135°C or colder (in liquid nitrogen). The cells are warrantied for twelve months from the date of receipt if stored at -135°C or colder. Storage of cells at -80°C for less than one month should maintain cell viability but is not covered by the warranty.

THAWING INSTRUCTIONS FOR CELL PRODUCTS

Refer to the Technical Support page at cgt.global/technical-support/ to access our online Thawing Protocols.

WARNING

This product contains human tissue or other biological material and MUST be handled at Biosafety Level 2 or higher. All biological products should be treated as potentially infectious or contaminated material, even if infectious disease screening reports are negative. Follow universal precautions and wear appropriate personal protective equipment.

PRODUCT WARRANTY

For our product warranty, please review our Terms and Conditions at cgt.global/terms-and-conditions/.

FOR RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.

CryoStor is a trademark of BioLife Solutions Inc

10032023-CGT-PIS-C-CB-CD8-V1 Page | 1