

1743 Creekside Drive, Suite 200 Folsom, CA 95630

530.626.7000 | www.cgt.global

email: info@cgt.global

# **Cryopreserved Product**

Peripheral Blood CD36+ Erythroid Progenitor Cells

Catalog# PB36001C 1.

1.0 million cells

PB36002.5C 2.5 million cells

#### PRODUCT DESCRIPTION

Human Peripheral Blood CD36+ Erythroid Progenitor Cells are derived from cultured peripheral blood CD34+ hematopoietic stem cells (HSCs).

Peripheral blood mononuclear cells are collected using the Spectra Optia® Apheresis System. CD34+ cells are selected from the mononuclear cell pool using immunomagnetic anti-CD34 microbeads. The isolated CD34+ cells are then cultured in StemSpan™ SFEM, a serum-free expansion media for hematopoietic cells, supplemented with SCF, EPO and IL-3, to generate CD36+ cells. Cultured cells are harvested after 10 days and checked for the expression of CD36.

Cells were obtained using Institutional Review Board (IRB) approved consent forms and protocols.

# **CRYOPRESERVATION**

Cryopreserved products allow for prolonged storage before use. Cell products contain 10% DMSO to minimize cell death during freezing. All cryopreserved products are stored in containers designed and tested for ultra-low temperatures at long time intervals. We normally ship cryopreserved items on dry ice, but can also use a cryoshipper at the customer's request

# SAMPLE COLLECTION AND PROCESSING

All samples are collected on-site at our Stem Cell Collection Center. Apheresis donors are transfused with ACD-A during the collection process. Samples are then quickly processed in our on-site laboratory to achieve maximum viability and quality. Cryopreserved cells are frozen at -1°C/minute in a -80°C freezer, and then transferred to liquid nitrogen.

Infectious disease testing for HIV, HBV, and HCV is performed on a sample of donor blood. Only samples with negative results within 90 days of collection are shipped unless approved by the customer. All testing is performed by a CLIA-certified lab.

# **FORMAT**

Stem/progenitor cells are frozen in StemSpan<sup>™</sup> + 10% DMSO. We can also use freezing media as specified by the customer.

# **STORAGE**

Cryopreserved cells should be maintained at -135°C or colder (in liquid nitrogen). The cells are warrantied for twelve months from the date of receipt if stored at -135°C or colder. Storage of cells at -80°C for less than one month should maintain cell viability but is not covered by the warranty.

#### THAWING INSTRUCTIONS FOR CELL PRODUCTS

Refer to the Technical Support page at cgt.global/technical-support/ to access our online Thawing Protocols.

# **WARNING**

This product contains human tissue or other biological material and MUST be handled at Biosafety Level 2 or higher. All biological products should be treated as potentially infectious or contaminated material, even if infectious disease screening reports are negative. Follow universal precautions and wear appropriate personal protective equipment.

# PRODUCT WARRANTY

For our product warranty, please review our Terms and Conditions at cgt.global/terms-and-conditions/.

FOR RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.

CryoStor is a trademark of BioLife Solutions Inc. Spectra Optia is a trademark of Terumo BCT, Inc.

10042023-CGT-PIS-C-PB-CD36-V1 Page | 1