



1743 Creekside Drive, Suite 200  
Folsom, CA 95630

530.626.7000 | [www.cgt.global](http://www.cgt.global)

email: [info@cgt.global](mailto:info@cgt.global)

## Cryopreserved Product

Peripheral Blood CD4+CD25+ Regulatory T Cells

Catalog#

PB425002C

2 million cells

### PRODUCT DESCRIPTION

Human Peripheral Blood CD4+CD25+ Regulatory T Cells are isolated from peripheral blood mononuclear cells in a multi-step process.

First, peripheral blood mononuclear cells are collected using the Spectra Optia® Apheresis System. Next, cells expressing CD8, CD14, CD16, CD19, CD20, CD36, CD56, CD66b, CD123, TCR $\gamma/\delta$  and CD235a are depleted from the mononuclear cell population using immunomagnetic particles leaving purified, untouched CD4+ T lymphocytes. CD25+ cells are then positively selected using immunomagnetic anti-CD25 particles from the untouched CD4+ T lymphocyte population, leaving highly enriched CD4+CD25+ regulatory T cells (T-regs).

Cells were obtained using Institutional Review Board (IRB) approved consent forms and protocols.

### CRYOPRESERVATION

Cryopreserved products allow for prolonged storage before use. Cell products contain 10% DMSO to minimize cell death during freezing. All cryopreserved products are stored in containers designed and tested for ultra-low temperatures at long time intervals. We normally ship cryopreserved items on dry ice, but can also use a cryoshipper at the customer's request.

### SAMPLE COLLECTION AND PROCESSING

All samples are collected on-site at our Stem Cell Collection Center. Apheresis donors are transfused with ACD-A during the collection process. Samples are then quickly processed in our on-site laboratory to achieve maximum viability and quality. Cryopreserved cells are frozen at -1°C/minute in a -80°C freezer, and then transferred to liquid nitrogen.

Infectious disease testing for HIV, HBV, and HCV is performed on a sample of donor blood. Only samples with negative results within 90 days of collection are shipped unless approved by the customer. All testing is performed by a CLIA-certified lab.

### FORMAT

Isolated cells are normally frozen in CryoStor™ CS10 (10% DMSO). We can also use freezing media as specified by the customer.

### STORAGE

Cryopreserved cells should be maintained at -135°C or colder (in liquid nitrogen). The cells are warranted for twelve months from the date of receipt if stored at -135°C or colder. Storage of cells at -80°C for less than one month should maintain cell viability but is not covered by the warranty.

### THAWING INSTRUCTIONS FOR CELL PRODUCTS

Refer to the Technical Support page at [cgt.global/technical-support/](http://cgt.global/technical-support/) to access our online Thawing Protocols.

### WARNING

This product contains human tissue or other biological material and MUST be handled at Biosafety Level 2 or higher. All biological products should be treated as potentially infectious or contaminated material, even if infectious disease screening reports are negative. Follow universal precautions and wear appropriate personal protective equipment.

### PRODUCT WARRANTY

For our product warranty, please review our Terms and Conditions at [cgt.global/terms-and-conditions/](http://cgt.global/terms-and-conditions/).

FOR RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.